

## CHROMATOGRAPHY REFRIGERATOR PLUS

### CRP Series Standard Features

- Temperature Range 0° - 10° C. Factory Preset at 4 ° C
- LR-5 Microprocessor Temperature Controller offers more precise temperature settings and maintenance over mechanical thermostats (see pg. 1)
  - Digital Display
  - Audible and Visual High & Low Temperature Alarm
  - Door Ajar Alarm
  - Battery Back-up for audible and visual alarms
  - Power Failure Alarm
  - Sensor Error Alarm
  - Remote Alarm Contacts
  - High & Low Non-volatile Temperature Memory with reset
  - Audible Alarm Mute Button
- 120 VAC Operation
- Accessible Exterior Power Switch
- One each interior GFI protected duplex electrical outlet
- Two each 2" access ports with spring loaded self closing, gasketed metal covers
- 4" Swivel Casters
- CFC-Free refrigerant, foam insulation and packaging
- Automatic Defrost Refrigerators with built in electric heater assisted condensate evaporator
- Interior lights are safety shielded and switch controlled
- Self-closing doors. Hold-open latches provided for sliding doors
- Energy efficient, double pane vacuum insulated glass doors with Solarban®
- Electrostatically coated wire shelves, adjustable in 1/2" increments
- Stainless Steel interior floor helps resist corrosion
- White acrylic enamel interior offers higher visibility
- Accu-test Certified. Double testing helps ensure accurate temperature performance

### CRP Accessories

- ( See Accessories Page 20)
- Additional Electrical Outlets
  - Additional Access Ports
  - Chromatography Mast
  - Additional Shelves
  - Half Shelves
  - Door Locks
  - 120 VAC Operated, Flush mounted, 4" Circular Chart Temperature Recorder
  - Solid Door Conversion protects light sensitive product



(CRP-47)  
SHOWN WITH OPTIONAL CHART RECORDER



(CRP-72)

CHROMATOGRAPHY PLUS														
Model	Cu. Ft.	Doors	Shelves	Usable Int. Dim.			Ext. Dim.			H.P.	Amps	BTU/hr	Weight	Ship Weight
				W"	D"	H"	W"	D"	H"					
CRP-26	26	1 swing	4	26	25	57	30	30	79"	1/3	8	2174	321	361
CRP-47	47	2 slide	8	50	24	56	54	30	79"	1/2	11	3217	478	518
CRP-49	49	2 swing	8	50	25	56	54	30	79"	1/2	11	3217	472	512
CRP-69	69	3 slide	12	75	24	56	79	30	80"	1/2	12	4314	664	714
CRP-72	72	3 swing	12	75	25	56	79	30	80"	1/2	12	4314	660	710

\* Exterior height does not include 5" for caster height

## CHROMATOGRAPHY REFRIGERATOR BASIC

### CRB Series Standard Features

- Temperature Range 0° - 10° C. Factory Preset at 4 ° C
- Mechanical Thermostat
- 120 VAC Operation
- Accessible Exterior Power Switch
- One each interior GFI protected duplex electrical outlet
- One each 2" access port with spring loaded self closing, gasketed metal cover
- Leveling Legs
- CFC-Free refrigerant, foam insulation and packaging
- Automatic Defrost Refrigerators with built in condensate evaporator
- Interior lights are safety shielded and switch controlled
- Self-closing doors. Hold-open latches provided for sliding doors
- Energy efficient, double pane vacuum insulated glass doors with Solarban®
- Electrostatically coated wire shelves, adjustable in 1/2" increments
- Stainless Steel interior floor helps resist corrosion
- White acrylic enamel interior offers higher visibility
- Accu-test Certified. Double testing helps ensure accurate temperature performance




(CRB-26)  
SHOWN WITH OPTIONAL SOLID DOOR CONVERSION



(CRB-47)

### CRB Accessories

- ( See Accessories Page 20)
- Additional Electrical Outlets
  - Additional Access Ports
  - Chromatography Mast
  - Additional Shelves
  - Half Shelves
  - Door Locks
  - Solid Door Conversion protects light sensitive product



CHROMATOGRAPHY BASIC SERIES														
Model	Cu. Ft.	Doors	Shelves	Usable Int. Dim.			Ext. Dim.			H.P.	Amps	BTU/hr	Weight	Ship Weight
				W"	D"	H"	W"	D"	H"					
CRB-23	23	1 swing	4	24	26	57	27	30	79	1/3	8	1640	312	352
CRB-26	26	1 swing	4	26	25	57	30	30	79	1/3	8	2174	321	361
CRB-33	33	2 slide	8	36	24	57	40	30	79	1/2	10	3217	404	444
CRB-35	35	2 swing	8	36	24	57	40	30	79	1/2	10	3217	413	453
CRB-47	47	2 slide	8	50	24	56	54	30	79	1/2	11	3217	478	518
CRB-49	49	2 swing	8	50	25	56	54	30	79	1/2	11	3217	472	512
CRB-69	69	3 slide	12	75	24	56	79	30	80	1/2	12	4314	664	714
CRB-72	72	3 swing	12	75	25	56	79	30	80	1/2	12	4314	660	710