

## DUAL TEMPERATURE UNITS FEATURING MANUAL DEFROST FREEZERS



(RFC-17A)

### RFC Series Standard Features

- Temperature Range:  
Refrigerator Compartment -  
Average Cycle Temperature of 4° C  
Freezer Compartment -  
Average Cycle Temperature Of -20° C or lower
- Mechanical Thermostat
- 120 VAC Operation
- Manual Defrost Freezer Compartments  
( Auto-Defrost on RFC-23A, RFC-49A & RFC-17A )
- Casters (Standard on RFC-49A, RFC-23A & RFC-13SS only)
- Door Locks (Standard on RFC-23A & RFC 49A only)
- CFC-Free refrigerant, foam insulation and packaging
- Two Digital Temperature displays on RFC-23A & RFC-49A
- Magnetic Door Gasket
- Accu-test Certified. Double testing helps ensure accurate temperature performance

### RFC Accessories

( See Accessories Page 20 )

- Caster Platform (available for RFC-34M only)
- Casters (available for RFC-18M)
- Battery Operated, Surface Mount, Digital Thermometer & Alarm (simple field installation required)
- Battery Operated Freestanding Internal 4" Circular Chart Temperature Recorder
- 120 VAC Operated, Freestanding Multi-probe temperature monitor & alarm; Great for combo units (simple field installation required)



(RFC-34M)

SHOWN WITH OPTIONAL CASTER PLATFORMS



(RFC-9M)



(RFC-9M)



(RFC-13SSM)



(RFC-49A)

### REFRIGERATOR/FREEZER COMBINATION UNITS

| Model     | Cu. Ft. |             | Defrost | Lights   | Shelves-Drawers | Usable Int. Dim. Ref/Fre |         |         | Ext. Dim. |         |         | H.P     | Amps | BTU/hr    | Ship |      |
|-----------|---------|-------------|---------|----------|-----------------|--------------------------|---------|---------|-----------|---------|---------|---------|------|-----------|------|------|
|           | Ref/Fre | Ref/Fre     |         |          |                 | Ref/Fre                  | Ref/Fre | Ref/Fre | Ref/Fre   | Ref/Fre | Ref/Fre |         |      |           | lbs. | lbs. |
| RFC-3M    | 2/1     | Auto/Manual | 2 swing | yes /no  | 2-1 / 0         | 16/14                    | 15/15   | 19/8    | 19        | 20      | 34      | 1/6     | 2    | 310       | 64   | 104  |
| RFC-9M    | 5.3/3.6 | Auto/Manual | 2 swing | yes /no  | 3-2 / 3-1       | 19/16                    | 18/18   | 29/27   | 22        | 24      | 66      | 1/4     | 7    | 850       | 127  | 167  |
| RFC-11M   | 8/3     | Auto/Manual | 2 swing | yes /no  | 5-1/1-0         | 21/19                    | 18/17   | 45/14   | 24        | 24      | 68      | 1/4     | 7    | 850       | 130  | 170  |
| RFC-13SSM | 6.5/6.5 | Auto/Manual | 2 swing | no /no   | 2-0 / 2-0       | 25/25                    | 22/23   | 24/24   | 27        | 29      | 77      | 1/5-1/3 | 5/4  | 770/1990  | 320  | 360  |
| RFC-18M   | 9/9     | Auto/Manual | 2 swing | yes /no  | 4-1 / 7         | 19/16                    | 18/16   | 51/50   | 44        | 24      | 57      | 1/6-1/6 | 7/7  | 850/850   | 220  | 260  |
| RFC-34M   | 17/17   | Auto/Manual | 2 swing | yes /yes | 4-0 / 4-0       | 25/25                    | 20/20   | 56/56   | 61        | 28      | 74      | 1/8-1/3 | 3/8  | 870/1988  | 350  | 390  |
| RFC-17A   | 13/4    | Auto/Auto   | 2 swing | yes /no  | 3-2/0/0         | 24/23                    | 19/19   | 39/15   | 28        | 32      | 65      | 1/4     | 7    | 850       | 180  | 220  |
| RFC-23A   | 11/8    | Auto/Auto   | 2 swing | yes /yes | 2-0 / 1-0       | 24/24                    | 26/26   | 30/20   | 27        | 30      | 79      | 1/5-1/3 | 12   | 770/1990  | 318  | 358  |
| RFC-49A   | 23/23   | Auto/Auto   | 2 swing | yes /yes | 3-0 / 3-0       | 24/24                    | 26/26   | 57/57   | 54        | 30      | 79      | 1/3-1/2 | 16   | 1990/4314 | 480  | 530  |