

GENERAL PURPOSE REFRIGERATORS & FREEZERS

GENERAL PURPOSE SOLID DOOR REFRIGERATORS

ARU Series Standard Features

- Temperature Range: 0° - 10° C. Factory preset at 4° C
- Mechanical Thermostat
- 120 VAC Operation
- CFC-Free refrigerant, foam insulation and packaging
- Interior Lights
- Full Length Door Handle
- Magnetic Door Gasket
- Accu-test Certified. Double testing helps ensure accurate temperature performance

ARU Accessories

(See Accessories Page 20)

- Caster Platform or Casters
- Digital Display & Alarm - Battery Operated
- Digital Display, Alarm & Remote Contacts, 115V operation
- Battery Operated, Digital Thermometer & Alarm (simple field installation required)
- Battery Operated Freestanding Internal 4" Circular Chart Temperature Recorder
- 120 VAC Freestanding Multi-probe temperature monitor & alarm; Great for combo units (simple field installation required)



(ARU-0904)



(ARU-1704)



MANUAL & AUTO DEFROST FREEZER BASIC

MFB, MCFB & AFB Series Standard Features

- Temperature Range: Average Cycle Temperature of -20° C or lower (MFB-2030 Average Cycle Temperature of -30° C)
- Mechanical Thermostat
- 120 VAC Operation
- Leveling Legs
- CFC-Free refrigerant, foam insulation and packaging
- Manual Defrost MFB & MCFB Units
- Auto-Defrost AFB Units
- Direct Freeze Shelves for faster temperature recovery (Excludes MCFB & AFB Units)
- Interior lights (Excludes MFB-0920)
- Door Locks
- Magnetic Door Gasket
- Accu-test Certified. Double testing helps ensure accurate temperature performance

MFB, MCFB & AFB Accessories

(See Accessories Page 20)

- Tower Racks
- Plastic Storage Boxes & Grids
- Cardboard Storage Boxes & Grids
- CryoGloves
- Battery Operated, Surface Mount Digital Thermometer & Alarm (simple field installation required)
- Battery Operated Freestanding internal 4" Circular Chart Temperature Recorder (for chest freezers)
- Casters or Caster Platforms



(MFB-0920)

(MCFB-2020)

(MFB-2020)

GENERAL PURPOSE REFRIGERATORS														
Model	Cu. Ft.	Doors	Shelves	Usable Int. Dim.			Ext. Dim.			H.P.	Amps	BTU/hr	Weight	Ship Weight
				W"	D"	H"	W"	D"	H"					
ARU-0904	9	1 swing	4/1 bin	19	18	51	57	22	24	1/6	7	850	94	134
ARU-1704	17	1 swing	4/3 bins	25	20	56	30	28	67	1/8	3	870	170	210
ARU-3104	31	1 swing	4/1 bin	32	26	65	35	32	72	1/3	5	2174	296	343

MANUAL/AUTO DEFROST BASIC																						
Model	Cu. Ft.	Defrost	Doors	Shelves	Usable Int. Dim.			Ext. Dim.			Space between fixed shelves/bins				H.P.	Amps	BTU/hr	lbs.	Ship lbs.			
					W"	D"	H"	W"	D"	H"	Top to Bottom											
MANUAL DEFROST BASIC																						
MFB-0920	9	Manual	1 swing	7 fixed	16	16	50	22	24	57	1.75	9	9.25	6.75	7	125	9.5	1/6	7	850	117	156
MFB-2020	20	Manual	1 swing	4 fixed	27	18	57	32	29	70	10.7	10.7	9.5	13.6	9.6			1/4	5	770	295	335
MFB-2030	20	Manual	1 swing	4 fixed	27	18	57	32	29	70	10.7	10.7	9.5	13.6	9.6			1/4	5	770	295	335
MFB-2920	29	Manual	1 swing	4 fixed/1 bin	32	26	65	35	32	77	11.8	13	12.75	12.75	9			1/3	5	2174	308	348
AUTO DEFROST BASIC																						
MCFB-1520	15	Manual	1 swing	2 baskets	42	24	28	48	30	35								1/4	5	770	161	206
MCFB-2020	20	Manual	1 swing	2 baskets	56	24	28	62	30	35								1/4	5	770	190	235
AFB-2020	20	Auto	1 swing	4 fixed	27	17	57	32	29	70	8.75	8.5	13.5	8.5	9.2			1/4	5	770	235	293
AFB-2920	29	Auto	1 swing	4 fixed/1 bin	32	26	66	35	32	72	12.75	12.8	12.75	12.75	12.75			1/2	9	2330	326	368