

LABORATORY REFRIGERATOR PLUS

LRP Series Standard Features

- Temperature Range 0° - 10° C. Factory Preset at 4° C
- LR-5 Microprocessor Temperature Controller offers more precise temperature settings and maintenance over mechanical thermostats (see pg. 1)
 - Digital Display
 - Audible and Visual High & Low Temperature Alarm
 - Door Ajar Alarm
 - Battery Back-up for audible and visual alarms
 - Power Failure Alarm
 - Sensor Error Alarm
 - Remote Alarm Contacts
 - High & Low Non-volatile Temperature Memory with reset
 - Audible Alarm Mute Button
- 120 VAC Operation
- Accessible exterior power switch
- 4" Swivel Casters
- CFC-Free refrigerant, foam insulation and packaging
- Automatic Defrost Refrigerators with built in electric heater assisted condensate evaporator
- Interior lights are safety shielded and switch controlled
- Self-closing doors. Hold-open latches provided for sliding doors
- Energy efficient, double pane vacuum insulated glass doors with Solarban®
- Electrostatically coated wire shelves adjustable in 1/2" increments
- Stainless Steel interior floor helps resist corrosion
- White acrylic enamel interior offers higher visibility
- Accu-test Certified. Double testing helps ensure accurate temperature performance



LRP Accessories

(See Accessories Page 20)

- Electrical Outlets
- Access Ports
- Additional Shelves
- Half Shelves
- Door Locks
- 120 VAC Operated, Flush mounted, 4" Circular Chart Temperature Recorder
- Solid Door Conversion protects light sensitive product

LAB REFRIGERATOR PLUS														
Model	Cu. Ft.	Doors	Shelves	Usable Int. Dim.			Ext. Dim.			H.P.	Amps	BTU/hr	Weight	Ship Weight
				W"	D"	H"	W"	D"	H"					
LRP-26	26	1 swing	4	26	25	57	30	30	79*	1/3	8	2174	321	361
LRP-47	47	2 slide	8	50	24	56	54	30	79*	1/2	11	3217	478	518
LRP-49	49	2 swing	8	50	25	56	54	30	79*	1/2	11	3217	472	512
LRP-69	69	3 slide	12	75	24	56	79	30	80*	1/2	12	4314	664	714
LRP-72	72	3 swing	12	75	25	56	79	30	80*	1/2	12	4314	660	710

* Exterior height does not include 5" for caster height

LABORATORY REFRIGERATOR BASIC

LRB Series Standard Features

- Temperature Range 0° - 10° C. Factory Preset at 4° C
- Mechanical Thermostat
- 120 VAC Operation
- Accessible Exterior Power Switch
- Leveling Legs
- CFC-Free refrigerant, foam insulation and packaging
- Automatic Defrost Refrigerators with built in condensate evaporator
- Interior lights are safety shielded and switch controlled
- Self-closing doors. Hold-open latches provided for sliding doors
- Energy efficient, double pane vacuum insulated glass doors with Solarban®
- Electrostatically coated wire shelves adjustable in 1/2" increments
- Stainless Steel interior floor helps resist corrosion
- White acrylic enamel interior offers higher visibility
- Accu-test Certified. Double testing helps ensure accurate temperature performance



LRB Accessories

(See Accessories Page 20)

- Electrical Outlets
- Access Ports
- Additional Shelves
- Half Shelves
- Door Locks
- Solid Door Conversion protects light sensitive product



LAB REFRIGERATOR BASIC														
Model	Cu. Ft.	Doors	Shelves	Usable Int. Dim.			Ext. Dim.			H.P.	Amps	BTU/hr	Weight	Ship Weight
				W"	D"	H"	W"	D"	H"					
LRB-12	12	1 swing	3	22	19	45	25	23	63	1/5	5	340	233	273
LRB-19	19	1 swing	4	24	21	57	27	25	75	1/3	5	1640	265	305
LRB-23	23	1 swing	4	24	26	57	27	30	79	1/3	8	1640	312	352
LRB-26	26	1 swing	4	26	25	57	30	30	79	1/3	8	2174	321	361
LRB-33	33	2 slide	8	36	24	57	40	30	79	1/2	10	3217	404	444
LRB-35	35	2 swing	8	36	24	57	40	30	79	1/2	10	3217	413	453
LRB-47	47	2 slide	8	50	24	56	54	30	79	1/2	11	3217	478	518
LRB-49	49	2 swing	8	50	25	56	54	30	79	1/2	11	3217	472	512
LRB-69	69	3 slide	12	75	24	56	79	30	80	1/2	12	4314	664	714
LRB-72	72	3 swing	12	75	25	56	79	30	80	1/2	12	4314	660	710