



# UNDERCOUNTER SERIES

Refrigerators & Freezers with front-breathing air-flow engineered specifically for built-in undercounter installations. This design allows for safe and efficient operation even when the unit is enclosed by laboratory and pharmacy casework.

## BUILT-IN REFRIGERATORS & FREEZERS



(UCBI-0620SSA)

### UCBI Series Standard Features

- Temperature Range:  
Refrigerators - Average Cycle Temperature of 4° C  
Freezers - Average Cycle Temperature of -20° C or Lower  
(UCBI - 0615 Average Cycle Temperature of -15° C)
- Mechanical Thermostat
- 120 VAC Operation
- Leveling Legs/Casters
- CFC-Free refrigerant, foam insulation and packaging
- Magnetic Door Gasket
- Accu-test Certified. Double testing helps ensure accurate temperature performance

### UCBI Accessories

( See Accessories Page 20 )


- Door Locks
- Battery Operated, Surface Mount, Digital Thermometer & Alarm (simple field installation required)
- 120 VAC Operated, Surface Mount, Digital Display, Alarm & Remote Contacts
- Battery Operated Freestanding Internal 4" Circular Chart Temperature Recorder
- 120 VAC Operated, Wall Mounted, 6" Circular Chart Temperature Recorder
- Casters or Caster Platforms
- Plastic Storage Boxes & Grids
- Cardboard Storage Boxes & Grids
- CryoGloves



(UCBI-0604)



(UCBI-0604SS)



UNDERCOUNTER REFRIGERATORS & FREEZERS																	
Model	Cu. Ft.	Temp.	Defrost	Doors	Base	Shelves/Bins	Usable Int. Dim.			Ext. Dim.			Ship				
							W"	D"	H"	W"	D"	H"	H.P.	Amps	BTU/hr	lbs.	lbs.
Front Breathing For Built-In Applications																	
UCBI-0604	6	+4°C	Auto	1 swing	Lev. Legs	2 removeable/0	20	17	26	24	24	34	1/5	3	310	100	140
UCBI-0615	6	-15°C	Manual	1 swing	Lev. Legs	3 removeable/0	20	17	26	24	24	34	1/5	3	310	100	140
UCBI-0604SS	6	+4°C	Auto	1 swing	Base	2 adjustable/0	25	23	24	27	30	34	1/5	3	310	215	255
UCBI0620SSA	6	-20°C	Auto	1 swing	Base	2 adjustable/0	25	23	24	27	30	34	1/3	5	2174	215	255

Compact Refrigerators & Freezers for undercounter "knee hole" applications where clearance is available to allow necessary airflow around external walls of unit.

## FREESTANDING REFRIGERATORS & FREEZERS

### UCFS Series Standard Features

- Temperature Range:  
Refrigerators - Average Cycle Temperature of 4° C  
Freezers - Average Cycle Temperature of -20° C or Lower
- Mechanical Thermostat
- 120 VAC Operation
- Leveling Legs
- CFC-Free refrigerant, foam insulation and packaging
- Magnetic Door Gasket
- Accu-test Certified. Double testing helps ensure accurate temperature performance

### UCFS Accessories

( See Accessories Page 20 )

- Door Locks
- Battery Operated Freestanding Internal 4" Circular Chart Temperature Recorder
- 120 VAC Operated, Surface Mount, Digital Display, Alarm & Remote Contacts
- 120 VAC Operated, Wall Mounted, 6" Circular Chart Temperature Recorder
- Casters or Caster Platforms
- Plastic Storage Boxes & Grids
- Cardboard Storage Boxes & Grids
- CryoGloves



(UCFS-0604SS)

(UCFS-0620SSM)



(UCFS-0320M)

(UCFS-RFC3M)

UNDERCOUNTER REFRIGERATORS & FREEZERS																	
Model	Cu. Ft.	Temp.	Defrost	Doors	Base	Shelves/Bins	Usable Int. Dim.			Ext. Dim.			H.P.	Amps	BTU/hr	lbs.	Ship lbs.
							W"	D"	H"	W"	D"	H"					
For Freestanding Or Kneehole Placement																	
UCFS-RFC3M	3			2 swing	Lev. Legs					19	20	34	1/6	2	310	64	104
Ref.	2	+4°C	Auto			1 adjustable/1	16	15	19								
Frz.	1	-20°C	Manual			0/0	14	15	8								
UCFS-0404	4	+4°C	Auto	1 swing	Lev. Legs	4 adjustable/0	18	17	28	21	22	34	1/6	2	310	70	110
UCFS-0320M	3	-20°C	Manual	1 swing	Lev. Legs	0/4	18	16	25	24	24	34	1/6	2	310	73	113
UCFS-0604SS	6	+4°C	Auto	1 swing	Casters	2 adjustable/0	25	22	24	27	29	33	1/5	5	770	135	185
UCFS-0620SSM	6	-20°C	Manual	1 swing	Casters	2 adjustable/0	25	23	24	27	29	33	1/3	4	1990	147	187